

# Cochlear™

## Nucleus® CP800 and CR100 Series

### Product Catalogue

*Hear now. And always*



**Cochlear™**


# Symbols



Note  
Important information or advice.  
Can avoid inconvenience.

---

## Colour index

	White
	Black
	Charcoal
	Brown
	Sand

This document shows all components and accessories available for the Cochlear™ Nucleus® CP800 Series Sound Processor and Cochlear Nucleus CR100 Series Remote Assistant.

To order items in this product catalogue please contact your clinic or nearest Cochlear office.

Please note that:

- Images shown are not actual size.
- Products are available as single items only, unless otherwise stated.
- Products shown in this catalogue might not be available in all countries and configurations or product specifications may vary.

# Contents

<b>Sound processor configurations</b>	<b>4</b>
CP810 configuration	4
CP800 LiteWear configuration	5
<b>CP810 Sound Processor - Kit</b>	<b>6</b>
<b>CR110 Remote Assistant - Kit</b>	<b>8</b>
<b>CP800 Series LiteWear - Kit</b>	<b>9</b>
<b>CP800 Series Rechargeable Battery System - Kit</b>	<b>10</b>
<b>CP800 Series Audio Accessories - Kit</b>	<b>11</b>
<b>Individual product components</b>	<b>12</b>
Processing unit	12
Remote assistant	12
Battery holder	13
Standard tamper resistant battery cover	13
Compact rechargeable battery module	13
Standard rechargeable battery module	14
Coils	14
Coil magnets	15
Coil spacer	16
Coil cables	16
Charging kits	17
<b>Accessories</b>	<b>18</b>
Sound processor covers	18
Snugfit™	19
Earhooks	19
LiteWear accessories	20
Remote assistant belt case	21
Remote assistant protective case	21
Remote assistant leather case	21
Freedom™ Accessory Adaptor	22
Monitor earphones	22
Audio accessories	23
FM cables	24
Zephyr by Dry & Store®	24
Dry-Brik® Dessicant	24
Everyday case	25
Everyday case drying capsules	25
Storage and accessories case	25
Microphone protector applicator	25
Zinc Air batteries	25

# Sound processor configurations

## CP810 configuration

- ① Coil
- ② Coil cable
- ③ Coil magnet
- ④ Processing unit
- ⑤ Earhook (standard)
- ⑥ Battery holder
- ⑦ Standard tamper resistant battery cover
- ⑧ Standard tamper resistant battery module
- ⑨ Standard rechargeable battery module
- ⑩ Compact rechargeable battery module



## CP800 Series LiteWear configuration

- ①① LiteWear cable
- ①② LiteWear case
- ①③ LiteWear fixing aid
- ①④ Snugfit™
- ①⑤ Tamper resistant earhook (small)



# CP810 Sound Processor - Kit

Kit number CN5-CP810

Please choose the colour



**Processing unit** \_\_\_\_\_

Standard earhook pre-attached



**Coil cable** (2 per kit) \_\_\_\_\_

Please choose the length

6 cm    8 cm    11 cm    28 cm



**Coil** \_\_\_\_\_



Magnet not included with coil (see page 14).



**Coil magnet** \_\_\_\_\_

Please choose the strength

½ M    1 M    2 M    3 M    4 M    5 M



**Battery holder** (white only) \_\_\_\_\_



**Standard tamper resistant battery cover** \_\_\_\_\_



**Earhook (standard)** (pack of 3) \_\_\_\_\_



**Zephyr by Dry & Store®** \_\_\_\_\_

Dry-Brik® included

Please choose the region

Europe

United Kingdom



**Monitor earphones** \_\_\_\_\_



To operate, monitor earphones must be connected to the Freedom™ Accessory Adaptor.



**Freedom™ Accessory Adaptor** \_\_\_\_\_



**Zinc Air batteries (pack of 6)** \_\_\_\_\_



**Everyday case** \_\_\_\_\_



**Everyday case drying capsules (4x)** \_\_\_\_\_



**Microphone protector applicator** \_\_\_\_\_

Includes 2 sets of replacement microphone protectors (4 mic protectors in total) mounted onto 2 applicators.

**Microphone protector removal tool** \_\_\_\_\_



**Tamper resistant locking tool** \_\_\_\_\_



# CR110 Remote Assistant - Kit

Kit number CN5-CR110

Remote assistant \_\_\_\_\_



Remote assistant leather case \_\_\_\_\_

Lanyard included, leather material



Remote assistant charging kit \_\_\_\_\_





# CP800 Series LiteWear - Kit

Kit number CN5-LW

Please choose the colour



**LiteWear cable 30 cm** \_\_\_\_\_

Please choose the length

30cm      50cm

**Tamper resistant earhook (small)** (pack of 6) \_\_\_\_\_

Fitting tool and 6 fitting pins included

**LiteWear case** (pack of 2) \_\_\_\_\_



Use with standard battery modules only.

**LiteWear fixing aids** \_\_\_\_\_

Pack includes hook & loop, safety pin and alligator clip

**Snugfit™** (pack of 2) \_\_\_\_\_

Please choose the size

Small      Medium      Large

**Koala soft toy** \_\_\_\_\_



# CP800 Series Rechargeable Battery System - Kit

Kit number CN5-RB

2 battery modules will be delivered as part of the kit

Please choose the colour



Compact rechargeable battery module \_\_\_\_\_



Standard rechargeable battery module \_\_\_\_\_



Sound processor battery charger (white only) \_\_\_\_\_



Sound processor charging kit (white only) \_\_\_\_\_



# CP800 Series Audio Accessories - Kit

Kit number CN5-AA

1 personal audio cable (unilateral or bilateral) will be delivered as part of the kit

**Personal audio cable** (unilateral – 3.5 mm) \_\_\_\_\_



**Bilateral personal audio cable** (3.5 mm) \_\_\_\_\_



**Mains isolation cable** (3 m, 3.5 mm) \_\_\_\_\_



Personal audio cable required to connect mains isolation cable to processing unit.



**Lapel microphone** \_\_\_\_\_



# Individual product components

## Processing unit

Standard earhook pre-attached

Z185556

Z179938

Z185555

Z185554

Z185553



## Remote assistant

Z159280



## Battery holder

Z198041



## Standard tamper resistant battery cover

Z185509

Z185471

Z185507

Z185505

Z185503



## Compact rechargeable battery module

Z202222

Z202219

Z202221

Z202220

Z202218



## Standard rechargeable battery module

Z196080

Z196076

Z196077

Z196073

Z196072



## Coils

Z209886

Z209880

Z209884

Z209882

Z206950



Magnet not included with coil.

## Coil magnets

½ M

1 M

2 M

3 M

Z218494

Z218500

Z218506

Z218511

Z218495

Z218501

Z218507

Z218512

Z218496

Z218503

Z218508

Z218513

Z218497

Z218504

Z218509

Z218514

Z218498

Z218505

Z218510

Z218515



4 M

5 M

Z218516

Z218521

Z218517

Z218522

Z218518

Z218523

Z218519

Z218524

Z218520

Z218525

## Coil spacer

Z208284

pack of 2



## Coil cables

6 cm

8 cm

11 cm

28 cm

Z207507

Z207512

Z207517

Z207523

Z207508

Z207513

Z207518

Z207524

Z207509

Z207514

Z207519

Z207525

Z207510

Z207515

Z207520

Z207526

Z207511

Z207516

Z207521

Z207527





## Sound processor battery charger

Z196071



Requires sound processor charging kit (Z206776).



## Sound processor charging kit

Z206776



## Remote assistant charging kit

Z218960



# Accessories

## Sound processor covers

pack of 2



**Freedom Pink**  
Z206751



**Freedom Blue**  
Z206752



**Freedom Silver**  
Z206753



**Floral**  
Z206754



**Leopard**  
Z206755



**Zebra**  
Z206756



**Swirl**  
Z206757



**Pretty**  
Z206758



**Pink Camo**  
Z206759



**Green Camo**  
Z206760



**Spots**  
Z206761



**Soccer**  
Z206763



Not for use with Snugfit™, compact rechargeable battery module or LiteWear cable.

## Snugfit™

**Small**

Z206770

**pack of 2**

**Medium**

Z206771

**Large**

Z206772



## Earhooks

**Tamper resistant (small)**

Fitting tool and 6 fitting pins included

Z206743

**pack of 6**



**Standard**

Z206742

**pack of 3**

**Large**

Z206745



## LiteWear cable

**30 cm**

Z206749

Z206747

Z209493

Z206748

Z209468

**50 cm**

Z246024

Z246020

Z246022

Z246021

Z246023



## LiteWear case

Z206750

**pack of 2**



Use with standard battery modules only.



## LiteWear fixing aids

Pack includes hook & loop,  
safety pin and alligator clip

Z207334



## Remote assistant belt case

Lanyard included, leather material

Z207975



## Remote assistant protective case

Lanyard included, rubber material

Z207981

Z217252

pack of 2



## Remote assistant leather case

Lanyard included, leather material

Z218393



## Freedom™ Accessory Adaptor

Z208295



The Freedom Accessory Adaptor is required to operate the monitor earphones and FM cables.

## Monitor earphones

Z208300



To operate, monitor earphones must be connected to the Freedom Accessory Adaptor.

## Personal audio cable

Z208289

**unilateral – 3.5 mm**



## Bilateral personal audio cable

Z208292

**3.5 mm**



## Mains isolation cable

Z208296

**3 m, 3.5 mm**



Personal audio cable required to connect mains isolation cable to processing unit.



## Lapel microphone

Z208299



## FM cables

FM cable type	60 cm	90 cm
FM2-E	Z77089	Z60369
FM10-E	Z77090	Z60370
FM12-E	Z77092	Z60372
FM14-E	Z77094	Z60374
FM15-E	Z77095	Z60375
FM18-E	Z77098	Z60377



To operate, FM cables must be connected to the Freedom Accessory Adaptor.

## Zephyr by Dry & Store®

Dry-Brik® included

### Region

Europe	Z210351
United Kingdom	Z210353



Global power adaptor supplied with Dry & Store®.



## Dry-Brik® Dessicant

Dry-Brik® Dessicant (for use with Zephyr Dry & Store)

FEZ00500

pack of 3





## Everyday case

Z206780



## Everyday case drying capsule

N94527



## Storage and accessories case

Z60821



## Microphone protector applicator

Includes 2 sets of replacement microphone protectors (4 mic protectors in total) mounted onto 2 applicators.

Z239979

## Zinc Air batteries

B13033

**pack of 6**



# Notes



# Cochlear™

**Cochlear Ltd** (ABN 96 002 618 073) 14 Mars Road, Lane Cove NSW 2006, Australia Tel: 61 2 9428 6555 Fax: 61 2 9428 6353  
**Cochlear Americas** 13059 E Peakview Avenue, Centennial CO 80111, USA Tel: 1 303 790 9010 Fax: 1 303 792 9025  
**Cochlear AG** European Headquarters, Peter Merian-Weg 4, CH - 4052 Basel, Switzerland Tel: 41 61 205 0404 Fax: 41 61 205 0405  
**Cochlear Deutschland GmbH & Co. KG** Karl-Wiechert-Allee 76A, D - 30265 Hannover, Germany Tel: 49 511 542 770 Fax: 49 511 542 7770  
**Cochlear Europe Ltd** 9 Weybridge Business Park, Addlestone Road, Addlestone, Surrey KT15 2UF, United Kingdom Tel: 44 1932 871 500 Fax: 44 1932 871 526  
**Nihon Cochlear Co Ltd** Ochanomizu-Motomachi Bldg, 2-3-7 Hongo, Bunkyo-Ku, Tokyo, 113-0033, Japan Tel: 81 3 3817 0241 Fax: 81 3 3817 0245  
**Cochlear (HK) Ltd** Unit 1810, Hopewell Centre, 183 Queens Road East, Wan Chai, Hong Kong SAR Tel: 852 2530 5773 Fax: 852 2530 5183  
**Cochlear (HK) Ltd** Beijing Representative Office Unit 2205 - 2207, Tower B, 91 Jianguo Road, Chaoyang District, Beijing 1000022  
P.R. China Tel: 8610 8599 9924 Fax: 8610 8599 9804  
**Cochlear Ltd** (Singapore Branch) 6 Sin Ming Road, #01-16 Sin Ming Plaza Tower 2, Singapore 575585 Tel: 65 6553 3814 Fax: 65 6451 4105  
**Cochlear Korea Ltd** 5F, Seong San BD, 1689-5, Seocho-dong, Seocho-gu, Seoul, Korea Tel: 82 2 533 4450 Fax: 82 2 533 8408  
**Cochlear Benelux NV** Schalienhoedreef 20i, B - 2800 Mechelen, Belgium Tel: 32 1579 5511 Fax: 32 1579 5500  
**Cochlear Italia SRL** Via Augusto Murri 45/L, 40137 Bologna, Italy Tel: 39 051 741 9811 Fax: 39 051 392 062  
**Cochlear France S. A. S.** Route de l'Orme aux Merisiers, Z.I. Les Algorithmes - Bât. Homère, 91190 Saint Aubin, France Tel: 33 811 111 993 Fax: 33 160 196 499  
**Cochlear Nordic AB** Konstruktionsvägen 14, SE - 43533 Mölnlycke, Sweden Tel: 46 31 335 14 61 Fax: 46 31 335 14 60  
**Cochlear Tibbi Cihazlar ve Sağlık Hizmetleri Ltd. Sti.** Cubuklu Mah. Bogazici Cad., Bogazici Plaza No: 6/1, Kavacik, TR - 34805 Beykoz-Istanbul, Turkey Tel: 90 216 538 5900 Fax: 90 216 538 5919  
**Cochlear Canada Inc** 2500-120 Adelaide Street West, Toronto, ON M5H 1T1, Canada Tel: 1 416 972 5082 Fax: 1 416 972 5083  
**www.cochlear.com**

Nucleus is a registered trademark of Cochlear Limited.  
Cochlear, the elliptical logo, Freedom and Snugfit are trademarks of Cochlear Limited.  
Dry & Store is a registered trademark of Ear Technology Corporation.  
Dry-Brik is a trademark of Ear Technology Corporation.

The statements made in this guide are believed to be true and correct in every detail as of the date of publication.  
However, specifications are subject to change without notice.

Products shown in this catalogue might not be available in all countries and configurations or product specifications may vary.  
Please contact your local Cochlear office for further information.

Hear now. And always

234564 ISS3 NOV09  
Printed in Switzerland